



Product Datasheet

AlaRS antibody (orb394685)

Description	Rabbit polyclonal antibody to AlaRS.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	alanyl-tRNA synthetase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	107 kDa
Uniprot ID	P49588
Dilution Range	WB:1:200-1:1000, IHC:1:20-1:200, IF/ICC:1:20-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using AlaRS antibody.



Immunohistochemical staining of paraffin-embedded tissue using AlaRS antibody.