



Product Datasheet

CD246 antibody (orb394687)

Description	Mouse monoclonal antibody to CD246.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	anaplastic lymphoma receptor tyrosine kinase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	176 kDa
Uniprot ID	Q9UM73
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant protein



Immunohistochemical staining of paraffin-embedded tissue