



Product Datasheet

alpha 1 Antichymotrypsin antibody (orb394688)

Description	Rabbit polyclonal antibody to alpha 1 Antichymotrypsin.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	P01011
Dilution Range	WB:1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human blood tis...



Immunohistochemical staining of paraffin...