



Product Datasheet

alpha 1 Antichymotrypsin antibody (orb394689)

Description	Mouse monoclonal antibody to alpha 1 Antichymotrypsin.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	48 kDa, observed 65 kDa
Uniprot ID	P01011
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human blood tis...



Immunohistochemical staining of paraffin...