


Product Datasheet

Alpha galactosidase A antibody (orb394692)

Description	Rabbit polyclonal antibody to Alpha galactosidase A.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	galactosidase, alpha
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	49 kDa
Uniprot ID	P06280
Dilution Range	WB:1:500-1:5000, IF/ICC:1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK-293 cells u...