


Product Datasheet

alpha N Catenin antibody (orb394694)

Description	Rabbit polyclonal antibody to alpha N Catenin.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	catenin(cadherin-associated protein), alpha 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	105 kDa
Uniprot ID	P26232
Dilution Range	WB:1:1000-1:10000, IP:1:200-1:2000, IF/ICC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y cells u...