



Product Datasheet

AMOT antibody (orb394701)

Description	Mouse monoclonal antibody to AMOT.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	angiomotin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	118 kDa
Uniprot ID	Q4VCS5
Dilution Range	WB:1:1000-1:10000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of GC-7901 cells u...



Immunohistochemical staining of paraffin...