



Product Datasheet

AMPK beta 2 antibody (orb394705)

Description	Rabbit polyclonal antibody to AMPK beta 2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	protein kinase, AMP-activated, beta 2 non-catalytic subunit
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	35 kDa
Uniprot ID	O43741
Dilution Range	WB:1:500-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...