



Product Datasheet

Amylase alpha antibody (orb394707)

Description	Rabbit polyclonal antibody to Amylase alpha.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	amylase, alpha 2B(pancreatic)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	58 kDa
Uniprot ID	P19961
Dilution Range	WB:1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of mouse skeletal ...



Immunohistochemical
staining of paraffin...