



Product Datasheet

ANO10 antibody (orb394713)

Description	Rabbit polyclonal antibody to ANO10.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	anoctamin 10
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	76 kDa
Uniprot ID	Q9NW15
Dilution Range	WB:1:500-1:5000, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...