


Product Datasheet

APOD antibody (orb394722)

Description	Mouse monoclonal antibody to APOD.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	apolipoprotein D
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	22-33 kDa
Uniprot ID	P05090
Dilution Range	WB: 1:5000-1:50000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human plasma ti...