


Product Datasheet

ATP13A1 antibody (orb394737)

Description	Rabbit polyclonal antibody to ATP13A1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	ATPase type 13A1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	109 kDa
Uniprot ID	Q9HD20
Dilution Range	WB: 1:1000-1:4000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver tis...