

Product Datasheet

ATP1A2-Specific antibody (orb394739)

Description Rabbit polyclonal antibody to ATP1A2-Specific.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen ATPase, Na⁺/K⁺ transporting, alpha 2(+) polypeptide

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 100 kDa


Uniprot ID [P50993](#)

Dilution Range WB: 1:500-1:1000; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...