



Product Datasheet

ATP5C1 antibody (orb394741)

Description	Mouse monoclonal antibody to ATP5C1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	ATP synthase, H ⁺ transporting, mitochondrial F1 complex, gamma polypeptide 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	33 kDa
Uniprot ID	P36542
Dilution Range	WB: 1:500-1:2000; IHC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical staining of paraffin...