



Product Datasheet

B23 antibody (orb394747)

Description	Mouse monoclonal antibody to B23.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	nucleophosmin(nucleolar phosphoprotein B23, numatrin)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	35-40 kDa
Uniprot ID	P06748
Dilution Range	WB: 1:5000-1:50000; IF: 1:50-1:200; IHC: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using B23 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using B23 antibody.