


Product Datasheet

BBS2 antibody (orb394753)

Description	Mouse monoclonal antibody to BBS2.
Species/Host	Mouse
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Bardet-Biedl syndrome 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	80 kDa
Uniprot ID	Q9BXC9
Dilution Range	WB: 1:500-1:2000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using BBS2 antibody.