



Product Datasheet

C11orf79 antibody (orb394782)

Description	Rabbit polyclonal antibody to C11orf79.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	chromosome 11 open reading frame 79
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Uniprot ID	Q9NX18
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells using C11orf79 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using C11orf79 antibody.