



Product Datasheet

C14orf138 antibody (orb394788)

Description	Rabbit polyclonal antibody to C14orf138.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	chromosome 14 open reading frame 138
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	45-50 kDa
Uniprot ID	Q9H867
Dilution Range	WB: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 cells using...



Immunohistochemical staining of paraffin...