



Product Datasheet

C14orf94 antibody (orb394789)

Description	Rabbit polyclonal antibody to C14orf94.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	chromosome 14 open reading frame 94
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	28kd
Uniprot ID	Q9H6D7
Dilution Range	WB: 1:500-1:3000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis
of PC-3 cells usin...



Immunohistochemical
staining of paraffin...