



Product Datasheet

C16orf84 antibody (orb394793)

Description	Rabbit polyclonal antibody to C16orf84.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	chromosome 16 open reading frame 84
Form/Appearance	liquid
Storage	PBS with 0.1% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	65-68 kDa
Uniprot ID	Q2VPK5
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using C16orf84 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using C16orf84 antibody.