


Product Datasheet

C17orf28 antibody (orb394794)

Description	Rabbit polyclonal antibody to C17orf28.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	chromosome 17 open reading frame 28
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa, 89 kDa
Uniprot ID	Q8IV36
Dilution Range	WB: 1:500-1:2000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of fetal human bra...