


Product Datasheet

C20orf46 antibody (orb394807)

Description	Rabbit polyclonal antibody to C20orf46.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	chromosome 20 open reading frame 46
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q9NUR3
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...