


Product Datasheet

CAPZA1 antibody (orb394837)

Description	Mouse monoclonal antibody to CAPZA1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	capping protein(actin filament) muscle Z-line, alpha 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	33kd
Uniprot ID	P52907
Dilution Range	WB: 1:1000-1:8000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using CAPZA1 antibody.