

Product Datasheet

Carboxypeptidase A3 antibody (orb394838)

Description Rabbit polyclonal antibody to Carboxypeptidase A3.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen carboxypeptidase A3(mast cell)

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 49kd


Uniprot ID **P15088**

Dilution Range WB: 1:500-1:2400; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse pancreas ...



Immunohistochemical staining of paraffin...