

Product Datasheet

Cardiac Troponin I antibody (orb394841)

Description Rabbit polyclonal antibody to Cardiac Troponin I.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen troponin I type 3 (cardiac)

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 26 kDa


Uniprot ID **P19429**

Dilution Range WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of rat heart tissue...



Immunohistochemical staining of paraffin...