


Product Datasheet

Casein Kinase 1 delta antibody (orb394843)

Description	Rabbit polyclonal antibody to Casein Kinase 1 delta.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	casein kinase 1, delta
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	47kd
Uniprot ID	P48730
Dilution Range	WB: 1:500-1:3000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...