



Product Datasheet

CD10,MME antibody (orb394854)

Description	Rabbit polyclonal antibody to CD10,MME.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	membrane metallo-endopeptidase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	100 kDa
Uniprot ID	P08473
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney ti...



MME antibody