



Product Datasheet

Integrin alpha 2b antibody (orb394861)

Description	Rabbit polyclonal antibody to Integrin alpha 2b.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	integrin, alpha 2b(platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	125 kDa
Uniprot ID	P08514
Dilution Range	WB: 1:1000-1:10000 IP:1:1000-1:10000 IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human blood tis...



Immunohistochemical staining of paraffin...