


Product Datasheet

CEACAM3-Specific antibody (orb394867)

Description	Rabbit polyclonal antibody to CEACAM3-Specific.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	carcinoembryonic antigen-related cell adhesion molecule 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P40198
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...