



Product Datasheet

CGNL1 antibody (orb394873)

Description	Rabbit polyclonal antibody to CGNL1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	cingulin-like 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	149kd
Uniprot ID	Q0VF96
Dilution Range	WB: 1:1000-1:4000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 cells u...



Immunohistochemical staining of paraffin...