

Product Datasheet

Choline kinase alpha-Specific antibody (orb394875)

Description Rabbit polyclonal antibody to Choline kinase alpha-Specific.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen choline kinase alpha

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 50 kDa


Uniprot ID **P35790**

Dilution Range WB: 1:200-1:2000 IP: 1:200-1:2000 IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of COLO 320 cells ...



Immunohistochemical staining of paraffin...