

Product Datasheet

Chymotrypsin like protease CTRL1 antibody (orb394879)

Description Rabbit polyclonal antibody to Chymotrypsin like protease CTRL1.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen chymotrypsin-like

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 28 kDa


Uniprot ID **P40313**

Dilution Range WB: 1:1000-1:10000 IHC: 1:20-1:200 IF/ICC: 1:10-1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse pancreas ...



Immunohistochemical staining of paraffin...