

Product Datasheet

CISD2-Specific antibody (orb394881)

Description Mouse monoclonal antibody to CISD2-Specific.

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen CDGSH iron sulfur domain 2

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG2a

Clonality Monoclonal


Purity ≥95% as determined by SDS-PAGE

MW 15 kDa


Uniprot ID [Q8N5K1](#)

Dilution Range WB: 1:2000-1:16000; IHC: 1:20-1:200; IF: 1:50-1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of fetal human bra...



Immunohistochemical staining of paraffin...