


Product Datasheet

CIT antibody (orb394882)

Description	Rabbit polyclonal antibody to CIT.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	citron(rho-interacting, serine/threonine kinase 21)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	230kd,177kd
Uniprot ID	O14578
Dilution Range	WB: 1:500-1:2000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells using CIT antibody.