


Product Datasheet

Claudin 18 antibody (orb394884)

Description	Mouse monoclonal antibody to Claudin 18.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	claudin 18
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	28-29 kDa
Uniprot ID	P56856
Dilution Range	IHC-P: 1:200-1200, IP: 1:200-1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...