



Product Datasheet

Claudin 4-specific antibody (orb394885)

Description	Rabbit polyclonal antibody to Claudin 4-specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	claudin 4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	22 kDa
Uniprot ID	O14493
Dilution Range	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse small int...



Immunohistochemical staining of paraffin...