



Product Datasheet

CNN2 antibody (orb394893)

Description	Mouse monoclonal antibody to CNN2.
Species/Host	Mouse
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	calponin 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	34-36kd
Uniprot ID	Q99439
Dilution Range	WB: 1:1000-1:4000; IHC: 1:500-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of pig stomach tis...



Immunohistochemical staining of paraffin...