



Product Datasheet

CPT1A antibody (orb394913)

Description	Mouse monoclonal antibody to CPT1A.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	carnitine palmitoyltransferase 1A(liver)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	86 kDa
Uniprot ID	P50416
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IP: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of L02 cells using...



Immunohistochemical staining of paraffin...