



Product Datasheet

CREST antibody (orb394922)

Description	Rabbit polyclonal antibody to CREST.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	synovial sarcoma translocation gene on chromosome 18-like 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	O75177
Dilution Range	WB: 1:500-1:2000; IP: 1:150-1:600; IHC: 1:20-1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using CREST antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using CREST antibody.