

Product Datasheet

Cryptochrome 2 antibody (orb394926)

Description Rabbit polyclonal antibody to Cryptochrome 2.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, IP, WB

Immunogen cryptochrome 2(photolyase-like)

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 65 kDa


Uniprot ID **Q49AN0**

Dilution Range WB: 1:500-1:1000; IP: 1:200-1:1000; IHC: 1:100-1:400; IF: 1:50-1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of Y79 cells using...



Immunohistochemical staining of paraffin...