



Product Datasheet

CSN1 antibody (orb394927)

Description	Rabbit polyclonal antibody to CSN1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	G protein pathway suppressor 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-70 kDa
Uniprot ID	Q13098
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using CSN1 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using CSN1 antibody.