



Product Datasheet

CTAGE1 antibody (orb394932)

Description	Mouse monoclonal antibody to CTAGE1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	cutaneous T-cell lymphoma-associated antigen 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	86-120 kDa
Uniprot ID	Q96RT6
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...