


Product Datasheet

CXCR3B-specific antibody (orb394936)

Description	Mouse monoclonal antibody to CXCR3B-specific.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	chemokine(C-X-C motif) receptor 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	40 kDa +43 kDa
Uniprot ID	P49682
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...