



Product Datasheet

CXorf61 antibody (orb394940)

Description	Rabbit polyclonal antibody to CXorf61.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	chromosome X open reading frame 61
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	13 kDa
Uniprot ID	Q5H943
Dilution Range	WB: 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...