


Product Datasheet

CYB5R3 antibody (orb394941)

Description	Mouse monoclonal antibody to CYB5R3.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	cytochrome b5 reductase 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	30-35 kDa
Uniprot ID	P00387
Dilution Range	WB: 1:1000-1:20000; IHC: 1:500-1:2000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...