



Product Datasheet

Cytoglobin antibody (orb394952)

Description	Mouse monoclonal antibody to Cytoglobin.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	cytoglobin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	21-29 kDa
Uniprot ID	Q8WWM9
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical staining of paraffin...