

Product Datasheet

Cytokeratin 19-specific antibody (orb394956)

Description Rabbit polyclonal antibody to Cytokeratin 19-specific.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen keratin 19

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 40 kDa


Uniprot ID [P08727](#)

Dilution Range WB: 1:500-1:1000; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of HepG2 cells using...



Immunohistochemical staining of paraffin...