

Product Datasheet

Cytokeratin 1-specific antibody (orb394957)

Description Rabbit polyclonal antibody to Cytokeratin 1-specific.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, IP, WB

Immunogen keratin 1

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 67 kDa


Uniprot ID [P04264](#)

Dilution Range WB:1:500-1:5000, IP:1:200-1:2000, IHC:1:20-1:200, IF/ICC:1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of mouse skin tiss...



Immunohistochemical staining of paraffin...