



Product Datasheet

Cytokeratin 7-specific antibody (orb394960)

Description	Rabbit polyclonal antibody to Cytokeratin 7-specific.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	keratin 7
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	54 kDa
Uniprot ID	P08729
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:2000; IHC: 1:20-1:500; IF: 1:20-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical staining of paraffin...