


Product Datasheet

DBC1 antibody (orb394965)

Description	Mouse monoclonal antibody to DBC1.
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	deleted in bladder cancer 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	88 kDa
Uniprot ID	O60477
Dilution Range	WB: 1:500-1:2000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of fetal human bra...